



# Immunohistochemistry (IHC)

DAB-Based Antibody Detection Protocol

## Protocol Overview

Immunohistochemistry (IHC) is a technique that uses enzyme-labeled antibodies to detect specific proteins or antigens in tissue sections. This protocol uses the ABC (Avidin-Biotin Complex) method with DAB (3,3'-Diaminobenzidine) chromogen for brown colorimetric detection.

### Duration

**Day 1:** ~4 hours (+ overnight incubation)

**Day 2:** ~4 hours

**Total:** 2 days

### Key Steps

1. Deparaffinization & antigen retrieval
2. Peroxide quench & blocking
3. Primary antibody (overnight)
4. Secondary antibody & ABC
5. DAB development & counterstain

### Important Notes:

- Prepare ABC solution at least 30 minutes before use
- Prepare DAB immediately before use - it's light-sensitive
- Monitor DAB development under microscope to prevent overstaining

- DAB is a potential carcinogen - handle with care and dispose properly

## Required Materials

### Reagents

- Xylene (for paraffin sections)
- 100%, 95%, 70%, 50%, 20% Ethanol
- ddH<sub>2</sub>O (distilled water)
- Citrate antigen retrieval buffer
  - Solution A: 0.1M citric acid
  - Solution B: 0.1M sodium citrate
  - Mix: 9 mL A + 41 mL B
- Filtered PBS
- PBS-T (PBS + 0.1% Triton X-100)
- 30% H<sub>2</sub>O<sub>2</sub> (hydrogen peroxide)
- 10% serum (matched to secondary antibody host)
- Primary antibody
- Secondary antibody (biotinylated)
- ABC kit (Avidin-Biotin Complex)
- DAB kit (3,3'-Diaminobenzidine)
- Hematoxylin (for counterstaining)
- 2% Acetic acid
- Bluing solution
- Histomount mounting medium

### Equipment

- Coplin jars
- Microwave
- Humidifying chamber
- PAP pen (hydrophobic barrier pen)
- Coverslips

- Timer
- 4°C refrigerator or cold room
- Light microscope (for monitoring DAB)
- Pipettes and tips
- Kimtech wipes

## Solution Recipes

### Citrate Buffer (pH 6.0):

Solution A: 19.21 g citric acid anhydrous in 1L H<sub>2</sub>O (or 21.01 g citric acid monohydrate)

Solution B: 29.41 g sodium citrate dihydrate in 1L H<sub>2</sub>O

Working solution: 9 mL A + 41 mL B

### 0.6% H<sub>2</sub>O<sub>2</sub> in PBS:

1 mL 30% H<sub>2</sub>O<sub>2</sub> + 49 mL filtered PBS

### PBS-T:

39 mL filtered PBS + 1 mL 10% Triton-X

### 10% Serum Block:

100 µL serum + 900 µL PBS-T

### 2% Serum (Antibody Dilution):

20 µL serum + 980 µL PBS-T

### DAB Solution:

2.5 mL H<sub>2</sub>O + 5 µL 30% H<sub>2</sub>O<sub>2</sub> + 1 drop buffer + 1 drop DAB reagent

⚠ Prepare immediately before use!

## Select Protocol Day:

Day 1 (4 hours + overnight)

Day 2 (4 hours)

## Day 1: Sample Preparation and Primary Antibody

Total time: ~4 hours + overnight incubation

Select Sample Type:

## 1 Deparaffinization (~30 minutes)

- ☐ Xylene - 10 minutes
- ☐ Xylene - 3 minutes
- ☐ Xylene - 3 minutes
- ☐ 100% EtOH - 2 minutes
- ☐ 100% EtOH - 2 minutes
- ☐ 95% EtOH - 1 minute
- ☐ 70% EtOH - 1 minute
- ☐ 50% EtOH - 1 minute
- ☐ 20% EtOH - 1 minute
- ☐ Transfer to ddH<sub>2</sub>O

## 2 Antigen Retrieval (~1 hour)

- ☐ Place slides in coplin jar with citrate buffer, sit 10 min at RT  
⚠ Have a second jar with buffer only to microwave alongside slides
- ☐ Microwave 2.5 min at 50% power (power level 5)

☐ Microwave 2.5 min at power level 4

☐ Microwave 2.5 min at power level 4

☐ Microwave 2.5 min at power level 4

☐ Remove lid and allow to cool AT LEAST 30 minutes

⚠ Critical: Do not rush cooling step

### 1 Antigen Retrieval (~15 minutes)

☐ Microwave citrate buffer 2.5 min at power level 5

☐ Pour heated buffer over slides

☐ Microwave slides 2.5 min at power level 4

☐ Allow to cool for 5 minutes

### 3 Peroxide Quench (~1 hour)

☐ Wash with filtered PBS 5 minutes

Pour citrate out of coplin jar, pour PBS in

☐ Prepare 0.6%  $\text{H}_2\text{O}_2$  in filtered PBS (1 mL 30%  $\text{H}_2\text{O}_2$  + 49 mL PBS)

💡  $\text{H}_2\text{O}_2$  is stored in 4°C fridge

☐ Treat slides with 0.6%  $\text{H}_2\text{O}_2$  for 30 minutes

Pour PBS out, pour  $\text{H}_2\text{O}_2$ /PBS in

- ☐ Wash 10 min with PBS-T 2× (in humidifying chamber)

PBS-T: 39 mL filtered PBS + 1 mL 10% Triton-X

#### 4 Blocking (1 hour)

- ☐ Prepare 10% serum in PBS-T (100 µL serum + 900 µL PBS-T)

💡 Use serum from same species as secondary antibody host (usually goat)

- ☐ Cover slides (approx. 300 µL per slide) with 10% serum for 1 hour

#### 5 Primary Antibody (~30 min + Overnight)

- ☐ Wash slides with PBS-T 5 min 3×

- ☐ Dilute primary antibody in 2% serum (20 µL serum + 980 µL PBS-T)

💡 Dilution varies by antibody - see Antibody Info tab

- ☐ Apply primary antibody to slides in humidifying chamber

- ☐ Incubate overnight at 4°C

⚠️ Ensure humidifying chamber is well-sealed

## Day 2: Detection, Development, and Mounting

Total time: ~4 hours

#### 1 Secondary Antibody (~1 hr 20 min)

- ☐ Wash slides with PBS-T 5 min 3×

- ☐ **Dilute secondary antibody in 2% serum**  
💡 Dilution varies by antibody - see Antibody Info tab

- ☐ **Overlay slides and incubate at room temp for 1 hour**

## 2 ABC Complex (~1 hr 20 min)

- ☐ **Prepare AB Solution AT LEAST 30 min before use**  
5 mL PBS-T + 1 drop A + 1 drop B, invert to mix, store in drawer

- ☐ **Wash slides with PBS-T 5 min 3×**

- ☐ **Overlay slides with AB solution, incubate at room temp 1 hour**

## 3 Color Development (DAB)

### ⚠️ CRITICAL STEP - DAB Safety

DAB is a potential carcinogen. Wear gloves, work in fume hood, and dispose properly as hazardous waste.

- ☐ **Wash slides with PBS-T 5 min 3×**

- ☐ **Prepare DAB IMMEDIATELY before use**  
2.5 mL H<sub>2</sub>O + 5 µL 30% H<sub>2</sub>O<sub>2</sub> + 1 drop buffer + 1 drop DAB reagent

- ☐ **Overlay slides with DAB, monitor under microscope**  
⚠️ Time and observe color development - typically 1-5 minutes

- ☐ **Place slides in H<sub>2</sub>O to stop reaction**

## 4 Counterstaining (~15 min)

☐ Place slides in hematoxylin for NO LONGER than 10 seconds

⚠ Over-staining will obscure DAB signal

☐ Wash slides with running H<sub>2</sub>O 5 minutes

☐ Dip slides in 2% acetic acid 10×

☐ Dip slides in ddH<sub>2</sub>O 10×

☐ Place slides in bluing solution 1 minute

☐ Wash slides with running H<sub>2</sub>O 3 minutes

## 5 Dehydration (~30 min)

☐ 20% EtOH - 1 minute

☐ 50% EtOH - 1 minute

☐ 70% EtOH - 1 minute

☐ 95% EtOH - 1 minute

☐ 100% EtOH - 1 minute

☐ 100% EtOH - 2 minutes

☐ 100% EtOH - 2 minutes

☐ Xylene - 5 minutes



☐ Xylene - 5 minutes

☐ Xylene - 5 minutes

## 6 Mounting and Coverslip

☐ Wipe xylene off slides with Kimtech wipes

⚠ DO NOT wipe tissue sample

☐ Add a couple drops of Histomount

☐ Place coverslip gently

☐ Use clean pipette tip to push out air bubbles

☐ Allow slides to dry overnight at room temperature

## Antibody Dilutions and Details

Below are specific antibody protocols used in our lab. Always optimize dilutions for your specific application.

### ATIC (5-Aminoimidazole-4-carboxamide ribonucleotide formyltransferase/IMP cyclohydrolase)

#### Storage & Location:

- Stored in -20°C freezer
- Antibody box A-D

#### Serum:

- Use goat serum for blocking

- Goat serum stored in -20°C freezer near TC

#### Dilutions:

- **Primary antibody:** 1:100 (10 µL per 1 mL solution)
- **Secondary antibody:** Goat anti-mouse, 1:200 (5 µL per 1 mL solution)
- Secondary stored in 4°C fridge, antibody box

## CC3 (Cleaved Caspase-3)

#### Storage & Location:

- Stored in -20°C freezer
- Antibody box A-D

#### Application:

- Apoptosis marker
- Commonly used for cell death studies

#### Dilutions:

Dilution to be optimized (protocol cut off in original document)

#### General Guidelines for New Antibodies:

- Always start with manufacturer's recommended dilution
- Run a dilution series (e.g., 1:50, 1:100, 1:200, 1:500) to optimize
- Include positive and negative control tissues
- Secondary antibody should be biotinylated for ABC method
- Match secondary antibody host to serum used for blocking
- Store antibody aliquots at -20°C to avoid freeze-thaw cycles

## Troubleshooting Guide

### High Background Staining

**Possible Causes:**

- Insufficient blocking
- Antibody concentration too high
- Inadequate washing
- Endogenous peroxidase not quenched
- DAB development too long

**Solutions:**

- Increase blocking time or serum concentration (try 20%)
- Increase antibody dilution (lower concentration)
- Extend wash times and increase number of washes
- Increase H<sub>2</sub>O<sub>2</sub> quenching time
- Monitor DAB more carefully under microscope

## No Signal or Weak Signal

**Possible Causes:**

- Antibody concentration too low
- Primary or secondary antibody not working
- Antigen not present or masked
- Insufficient DAB development
- ABC solution not prepared correctly

**Solutions:**

- Decrease antibody dilution (higher concentration)
- Test antibodies with positive control tissue
- Optimize antigen retrieval (try different methods)
- Extend DAB development time
- Ensure ABC solution was prepared 30+ min before use

## Excessive Counterstain

**Possible Causes:**

- Hematoxylin incubation too long

- Hematoxylin solution too concentrated

**Solutions:**

- Reduce hematoxylin time to 5-10 seconds
- Use fresh hematoxylin or dilute if too strong
- Can differentiate in acid alcohol if caught early

## DAB Precipitate on Slides

**Possible Causes:**

- DAB solution not mixed well
- DAB solution prepared too far in advance
- Dust or debris on slides

**Solutions:**

- Mix DAB solution thoroughly
- Always prepare DAB immediately before use
- Filter DAB solution if necessary
- Ensure slides are clean and dust-free

## Tissue Falls Off Slides

**Possible Causes:**

- Slides not positively charged
- Microwave power too high
- Excessive agitation during washes

**Solutions:**

- Use positively charged slides only
- Reduce microwave power during antigen retrieval
- Be gentle during washing steps
- Bake slides at 60°C for 30-60 min before staining



## Safety Information

- **DAB is a potential carcinogen** - always wear gloves, work in fume hood, dispose as hazardous waste
- Always work with xylene in a fume hood - fumes are toxic
- Wear appropriate PPE (lab coat, gloves, safety glasses)
- Handle antibodies with care - some may contain sodium azide preservative
- $H_2O_2$  can cause burns - handle with care
- Be cautious with hot liquids during antigen retrieval
- Follow proper storage conditions for antibodies (typically 4°C or -20°C)

[↑ Back to Top](#)

## About This Protocol

Immunohistochemistry (IHC) staining protocol using ABC method with DAB chromogen for detecting specific proteins in tissue sections.

## Quick Links

[Protocol Overview](#)

[Antibody Info](#)

[Troubleshooting](#)

## Contact

Román A Cáceres

University of Cincinnati:  
[Cacerera@mail.uc.edu](mailto:Cacerera@mail.uc.edu)

Emory University:  
[Racacer@emory.edu](mailto:Racacer@emory.edu)

